

TOWER RECIRCULATION

- Improved Air Quality = Higher Fertilization Rates
- 100 % Air recirculation in IVF lab
- Carbon filter will remove VOC & CAC free
Class 10,000 environment
- HEPA filters will remove 0.3 micron particles
- Pre Filter will remove 10 and 5 micron particles
- Blowers(2 Nos) with Noise level < 52dBA
- Standby & running mode.
- Magnehelic gauge will represent filter life time
- Easy movable with wheels
- Covers upto 200 Sq.Ft area
- Optional : Germicidal UVC



MOBILE AIR



- Movable to required area
- Activated Carbon impregnated with alumina & KOH removes molecular contaminants
- HEPA filtration removes particulate contaminants
- 4 speed option for variable room size.
- 51 dBA noise level so low in the industry
- Maximum Air Flow: 300 Cubic Feet/minute
- Air Exchanges/Hour in 150 sq. ft.: 9.53
- Effective Coverage: 150 sq. ft.
- 100% recirculation

A Product From Your Friend

Esquire IVF Equipments Pvt Ltd