

## DRY INCUBATOR



- Body made of GI with powder coat finish
- Chamber made of SS304
- Omron PID controller for better accuracy
- Two shelves for storage your disposable and semen samples
- Temperature range from 49.9°C
- Temperature accuracy :  $\pm 0.3^{\circ}\text{C}$
- Compact design
  - Chamber size 420 x 420 x 400mm
- Visible polycarbonate front door to view disposable and samples

*A Product From Your Friend*

**Esquire IVF Equipments Pvt Ltd**

5A, Sri Mahalakshmi Co-Operative Nagar, Vazhuthalampedu, Kundrathur, Chennai 600 069  
www.esquires.in | admin@esquires.in