



Workstation Incubator is a compactable incubator which has 4 dish chamber with premixed gases of CO_2 , N_2 and O_2 .

The Workstation Incubator provides a good solution for short time culture management during oocyte selection & denudation

The Workstation Incubator features:

- Temperature range ambient to 40°C
- Temperature accuracy $\pm 0.2^\circ\text{C}$
- Analog flow meter for adjusting flow rate of pre mixed gas.
- Optimal Digital flow meter
- Flushing and purge function for fast recovery in the chamber.
- Chamber with space for 2 culture dishes.
- Alarm for gas flow interruption.
- Temperature can be easily controlled
- Dish Plate is removable & autoclavable.
- Can be used as warming chamber too.

A Product From Your Friend

Esquire IVF Equipments Pvt Ltd