

STAR LAF 2 FEET



- Made of GI with powder coat finish
- Working Table S.S. 304
- Working area : 650 L x 600 W x 600 H in mm .
- Bottom rack for keeping accessories
- Vertical down flow lamhinar system
- HEPA filters with rating of 99.99% efficacy 0.3 μ particulate
- EBM german blower for continues operation.
- Illumination light & germicidal UV Incorporated
- Additional power point inside the working area

A Product From Your Friend

Esquire IVF Equipments Pvt Ltd

5A, Sri Mahalakshmi Co-Operative Nagar, Vazhuthalampedu, Kundrathur, Chennai 600 069
www.esquires.in | admin@esquires.in