


- Model: MS HB / SS HB
- Made of GI with powder coat finish / Made of stainless steel body
- Ormon PID controller
- Temperature accuracy $\pm 0.2^{\circ}\text{C}$
- Temperature range ambient to 49.9°C
- 3 No's of anodized aluminium blocks
- Each blocks can hold 6 No's of ASP tubes
- Optional 30° angle block
- Easy to clean block



A Product From Your Friend

Esquire IVF Equipments Pvt Ltd

#5A, Sri Mahalakshmi Co-Operative Nagar, Vazhuthalampedu, Kundrathur, Chennai 600 069
www.esquires.in | admin@esquires.in