

# INCUBATOR WORKING CHAMBER

## Product Highlights

IncuWork *ic*

PULSE SCIENCE

IncuWork is an IVF Working chamber, totally closed embryo culture system, which substantially enhances the working environment for IVF procedures.

### SPECIFICATIONS

#### External workstation dimensions

- Height (excluding stand): 725 mm
- Width (side to side): 1180 mm
- Depth: 690 mm

#### Internal workstation dimensions

- Height: 590 mm
- Width: 1174 mm
- Depth: 570 mm

#### Hatch box interlock dimension


- Height: 220 mm
- Width: 220 mm
- Depth: 320 mm

#### Atmosphere Control

- Temperature: ambient 5°C- 45°C
- O2 Control: 0 - 10%
- CO2 Control: 0 - 10%

#### Alarms

- Temperature: visual and audible
- Gas level: visual and audible




TOUCH-SCREEN DISPLAY

The IncuWork *ic* is designed for efficiency and comfort allowing the operation to reduce stress and fatigue. Easy and comfortable performance.

This product can reduce embryo stress, improve embryo implantation rates, improve egg quality and increase the clinical pregnancy rates.

*Ergonomic design*

*Stable O2, CO2, temperature and humidity*



ESQUIRE

*A Product From Your Friend*

Esquire IVF Equipments Pvt Ltd

#5A, Sri Mahalakshmi Co-Operative Nagar, Vazhuthalampedu, Kundrathur, Chennai 600 069  
www.esquires.in | admin@esquires.in