

Synvivo

Tabletop Multi-Chamber Incubator



Features

1. Excellent efficiency with 8 chambers
2. Heating lid/bottom of each chamber offers the best temperature environment
3. Individual gas management ensures minimum gas waste
3. User-friendly Graphical User Interface
4. All Log Data will be maintained and ready for download
5. Comfortable LCD Touch Monitor
6. Temperature Calibration, Alarm List, VOC Filter Management



ESQUIRE

A Product From Your Friend

Esquire IVF Equipments Pvt Ltd

#5A, Sri Mahalakshmi Co-Operative Nagar, Vazhuthalampedu, Kundrathur, Chennai 600 069
www.esquires.in | admin@esquires.in


TABLETOP INCUBATOR

Synvivo

Eight Chambers are all independently controlled



- ✓ Each chamber's incubation on/off and setting
- ✓ Each chamber's lapsed duration of incubating
- ✓ Each chamber's set and ongoing temperatures
- ✓ Icon for setting configuration
- ✓ Icon for temperature control
- ✓ List of triggered alarms
- ✓ GAS supply pressure
- ✓ VOC Filter expiration



ESQUIRE

A Product From Your Friend

Esquire IVF Equipments Pvt Ltd

#5A, Sri Mahalakshmi Co-Operative Nagar, Vazhuthalampedu, Kundrathur, Chennai 600 069
www.esquires.in | admin@esquires.in

TABLETOP INCUBATOR